



At least 48 hours before you travel, tell your airline operator Wizzybug's battery type [Lithium Ion], the energy rating [256 Wh] and how to prevent Wizzybug from switching on accidentally, as described here.

Specification	
Manufacturer: Designability Charity Ltd.	Manufacturer contact number: +44 (0)1225 824103
Product type: Electrically powered child's wheelchair	Battery type: 1 x Lithium Iron Phosphate [LiFePO ₄ , Lithium Ion]
Product classification/approval: Class I Medical Device	Battery safety: UN38.3 approved, internal short-circuit protection
Product dimensions: Approx. 650mm [L] x 550mm [W] x 850mm [H]	Battery energy rating: 256 Wh [Voltage 25.6 V, Capacity 10 Ah]
Product weight: 30.8 kg [incl. battery], 28.0 kg [excl. battery]	Battery weight: 2.8 kg

Owner Information	
Name of child this Wizzybug belongs to:	Name[s] of adult passenger[s]:
Outbound journey [airport] to [airport]:	Outbound flight number and date of travel:
Return journey [airport] to [airport]:	Return flight number and date of travel:
Contact name and telephone number:	Notes:

Turn over for 'How to Disconnect the Battery' →



How to Disconnect the Battery	Consult Instructions for Use for full guidance.
Wizzybug must be prevented from operating accidentally by disconnecting the battery.	
<p>1. Remove the seat by turning the white dial, then squeezing the knob and lifting by the harness.</p> 	<p>2. Open the battery cover by pinching the black tab and white button at the front of the cover.</p> 
<p>3. Disconnect the battery by sliding back the silver slider [this is part of the connector] and turning the connector anti-clockwise.</p> 	<p>4. As advised by the airline; close the cover and replace the seat to transport Wizzybug whole, or remove the battery to carry on board.</p> 

How to Move and Store Wizzybug	
<p>Lifting of the Wizzybug is a 2-person operation. Never lift by the headrest or armrests. Use the front wheel arches and black rear metal beam.</p> 	<p>To freewheel, slide the metal tab away from the green dot on both wheels. To stop freewheeling, slide the tab to the green position and push.</p> 

Turn over for 'Specification and Owner Information' →



Wizzybug Passport

The **Wizzybug Passport** on the next page is for a Wizzybug with a LiFePO4 battery. If you have an SLA or NiMH battery, please download and fill in a copy of **FRM-013 Wizzybug Travel**

Transport information for the Wizzybug

Wizzybug General Information: Weight: 31.3 kg Length: 650mm Width: 560mm Battery Type:.....

This **Wizzybug powered wheelchair** belongs to

Travelling with adults Phone number

Flight Information: From: Airport To: Airport Flight no: Date:/...../.....

Lifting



See above picture for lifting Wizzybug. It is a **2 person** lift. Never lift by the headrest or armrests. Instead use the front wheel arches and black back metal beam.

© Designability 2016
Registered Charity No. 256335

Wizzybug Travel Document

FRM-013 1.0 13th September 2016

Moving and Handling

To **freewheel**, slide the metal tab to the other end of the recess on both front wheels.

To **re-engage**, slide the metal tab back to the original position, then push Wizzybug forwards to re-engage the motor.

